



**Biocompatibility**  
Class II



**Biocompatibility**  
Class II



**Biocompatibility**  
Class I



**Library Updates**  
Via WiFi



**Library Updates**  
Via SD card

-



**SOL Integration**

-

-



**Door Opening Mechanism**  
One-Touch



**Door Opening Mechanism**  
Manual



**Door Opening Mechanism**  
Manual



**Interface**  
Touchscreen



**Interface**  
Non Interactive



**Interface**  
Non Interactive



**Glass Panel**  
Dynamic



**Glass Panel**  
Static



**Glass Panel**  
Static



**Chamber Dimensions**  
20(L) x 15(W) x 6(H) cm  
7.9 (L) x 5.9 (W) x 2.4 (H) in



**Chamber Dimensions**  
14(L) x 14(W) x 8(H) cm  
5.5 (L) x 5.5 (W) x 3.1 (H) in



**Chamber Dimensions**  
15 (L) x 17 (W) x 7 (H) cm  
5.9 (L) x 6.7 (W) x 2.8 (H) in



**Full Arches**  
8



**Full Arches**  
4



**Full Arches**  
4



**LED Lights**  
65



**LED Lights**  
31



**LED Lights**  
12