

CURIE^{Plus}

Professional UV Curing Made Effortless



High Capacity Curing

- Fits up to 8 full arches.



Validated with Top Resins

- Over 150+ validated resins.
- Class I & II validated industry-leading resins.



Compatible with Nitrogen Sources

- Built for nitrogen tanks and generators.
- Up to a 25% boost in mechanical strength with nitrogen curing.



No Flipping Prints Required

- Mix-matrix 365, 385, and 405 nm LEDs provide 360° curing coverage.

Item Number

CURIE Plus RS110010162

A New Standard in Simplicity

- One-touch open/close door and sliding glass platform.
- Seamless integration with your SOL and SOL Plus.



Validated Materials Available



View a complete list of validated materials.



Dental Model
Curing Time:
10 min



Temporary Crown
Curing Time:
5 min



Splint
Curing Time:
5 min



Denture Base
Curing Time:
10 min



Denture Teeth
Curing Time:
5 min



Surgical Guide
Curing Time:
6 min



IBT
Curing Time:
10 min



Castable Crown
Curing Time:
5 min

Technical Specifications



Flash Curing
10 flashes per second



UV Wavelengths
365, 385, 405 nm



LED Lights
65



Chamber Dimensions
7.9 (L) x 5.9 (W) x 2.4 (H) in
20 (W) x 15 (D) x 6 (H) cm



Dimensions
14.4 (W) x 10.2 (D) x 9.2 (H) in
36.6 (W) x 26 (D) x 23.4 (H) cm



Weight
26.9 lbs
12.2 kg



Power Consumption
500 W max.



Power Input
100-240VAC, 50-60Hz



Library Updates
Automatic online

Package Contents

1 x CURIE Plus UV Oven
1 x N2 Cover
1 x Glass Panel

1 x Power Cable
1 x Ethernet Cable
1 x Wi-Fi Dongle



Learn More



Endless Nitrogen, Enhanced Results



Item Number

NFINITY RS110010163



Up to 25% Increased Flexural Strength

- Enhanced color stability and surface smoothness.



Eliminate Tank Hassles

- Eliminates the need for ordering, storing, and replacing nitrogen tanks.



Space-Efficient

- Compact design integrates seamlessly into your existing workspace.



Never Run Out of Nitrogen

- Generate unlimited 95% pure nitrogen on demand.

Technical Specifications



Dimensions

11.6 (W) x 9.4 (D) x 14.2 (H) in
29.5 (W) x 24 (D) x 36 (H) cm



Power Input

100-220VAC, 50-60Hz



Nitrogen Flow Rate

Up to 1.6 L/min



Weight

37.9 lbs
17.2 kg



Compatibility

Nitrogen parameters
preset in CURIE Plus



Nitrogen Purity

≥ 95%

Package Contents

1 x NFINITY Nitrogen Generator

1 x Power Cable

1 x Nitrogen Tube

1 x Drainage Tube

2 x N2 Covers

2 x Lock Screws

1 x Magnetic Holder

3 x Tube Holders



Learn More