

CURIE

Versatility that Delivers



Efficient Curing Space

- Perfect footprint for dental offices and labs.



Future-proof

- Easy resin library update mechanism with SD card reader.



Validated Workflows

- Over 150 preloaded resin profiles available at the touch of a button.



Premium Performance

- 360-degree light distribution covering 365, 385 and 405 nm achieving superior curing results.



CURIE Technical Specifications



Curing Time Cycles

1 sec. to 20 min.



Dimensions

7.9 (W) x 8.7 (D) x 6.3 (H) in



Weight

6.6 lbs



Chamber Dimensions

5.5 (W) x 5.5 (D) x 3.1 (H) in



UV Wavelength

365, 385, 405 nm



Technology

UV LED Array



Power Consumption



75 W



Power Input

100-240VAC, 50-60Hz

Item Number
CURIE - KT200200101



PACKAGE CONTENTS

- 1 x Glass Panel Box
- 1 x Power Cable
- 1 x Ethernet Cable
- 1 x Wi-Fi Dongle



Scan the QR Code to
Schedule a Virtual Demo Today!

ackuretta.com
northamerica.office@ackuretta.com