



	CURIE	Otoflash	Meccatronicore	Uvitron Intelliray
Model	10 min.	2000 flashes per side	6 min.	2 min per side
Element Model	10 min.	2000 flashes per side	6 min.	2 min per side
Crown	10 min.	2000 flashes per side	6 min.	2 min per side
Splint	10 min.	2000 flashes per side	6 min.	2 min per side
Guide	10 min.	2000 flashes per side	6 min.	2 min per side
Tray	7 min.	2000 flashes per side	6 min.	2 min per side
ProCast	5 min.	2000 flashes per side	6 min.	2 min per side
Denture	5 min.	2000 flashes per side	6 min.	2 min per side
IBT	10 min.	2000 flashes per side	6 min.	2 min per side
Gingiva	/	2000 flashes per side	/	/