

SOL^{SE} | The Smart Choice for Speed and Precision



A New Benchmark in Production Speed

Print up to 8 arches nested vertically

Validated, Versatile Solution

Over 200 validated dental resins

Reliable and Efficient

Multiple build platform sizes available

2-year Priority Support

Direct Support with Concierge Service

Item Number

LD001 - SOL SE 3D Printer

LD002 - with Concierge Service*

*Available in the USA, Canada and Australia

Technical Specifications

Print Area

5.1 (W) x 3.1 (D) x 4.7 (H) in
12.9 (W) x 8 (D) x 12 (H) cm

Dimensions

11.5 (W) x 10.7 (D) x 16.4 (H) in
29.2 (W) x 27.1 (D) x 41.6 (H) cm

3D Print Technology

Monochrome LCD

Power Input

100-240VAC, 50-60Hz

UV Wavelength

385-405 nm

Layer Thickness

30-150 µm

Weight

40 lbs / 18 kg

XY Resolution

49 µm



Validated Materials Available

ARGEN

BEGO

DeltaMed
TURNING IDEAS INTO MATERIALS

dentona
NEW WAY DENTAL TECHNOLOGY

keystone
industries

NextDent

PacDent
Passion for Excellence

VOCO

Learn more about
SOL SE



Dental Model
Printing Time:
L 17 min



Permanent Crown
Printing Time:
L 10 min



Splint
Printing Time:
L 20 min



Denture Base
Printing Time:
L 25 min



Denture Teeth
Printing Time:
L 10 min



Surgical Guide
Printing Time:
L 20 min



IBT
Printing Time:
L 20 min



Castable Crown
Printing Time:
L 20 min

Package Contents

SOL SE 3D Printer

L Build Platform

Vat & Vat Box

Power Cable

Ethernet Cable

Wi-Fi Dongle

Hex Keys (2)

Trays (2)

Mallet

Spatula

Clipper

Scraper

Concierge Service Box*

*Only available with SOL SE
Concierge Service customer

Available Accessories



L Build Platform
Item Number
KT210010106



M Build Platform
Item Number
KT210050106



S Build Platform
Item Number
KT210040106



Vat Assembly with Vat Box
Item Number
RS110010166



Vat Box (2pcs)
Item Number
RS110010165



Print Screen (2pcs)
Item Number
KT110010106