

The Ackuretta **SOL** 3D Printing Solution

Premium Performance meets Affordability




Precision Chairside Production

- Print up to 8 full arches vertically & 3 full arches flat
- 49 µm XY pixel size creates smooth & fine prints reliably
- Variable layer thickness from 30 to 150 µm



Validated Workflows

- 150+ validated resin profiles available for use with slicing software and 3D printer, and curing unit
- Class I & Class II biocompatible workflows approved for a wide variety of FDA and Health Canada-certified resins



Reliable and Efficient

- 3 Build platform sizes available (Small, Medium and Large) to match your production needs
- Small build platform is ideal for fast crown, bridge, inlay and onlay prints
- Medium build platform prints fast full arch prints ideal for full arch implants, surgical guides, splints & more
- Large build platform is ideal for printing multiple cases efficiently



Continued Support & Training

- Concierge Service including remote onboarding, training & support throughout 2-year 3D printer warranty
- 2-year SOL warranty



Item Numbers

SOL with Concierge Service - RS110010240

CLEANI with Dual Tanks - RS110010246

CURIE Plus with N2 Cover - RS110010162



Scan the QR Code to
Schedule a Virtual Demo Today!

ackuretta.com
northamerica.office@ackuretta.com



Validated Materials Available in North America



For a complete list of all validated materials for Ackuretta's SOL 3D Printer, scan this QR Code.



SOL Technical Specifications



3D Print Technology
LCD



Dimensions
11.5 (W) x 10.7 (D) x 16.4 (H) in



Weight
40 lbs



Print Area
5.1 (W) x 3.1 (D) x 4.7 (H) in



XY Resolution
49 μ m




Layer Thickness
30-150 μ m




UV Wavelength
385-405 nm



Power Input
100-240VAC, 50-60Hz



SOL PACKAGE CONTENTS

- 1x SOL 3D Printer
- 1x  Build Platform
- 1x Ethernet Cable
- 1x Power Cable
- 1x WiFi-Dongle
- 2x Hex Key
- 1x Vat & Vat Storage Box
- 1x Starter Tool Kit
- 1x Concierge Service Box



Scan the QR Code to
Schedule a Virtual Demo Today!